




Paging Interface Adapter Enterprise System *PIAES*

*comPPage's **PIA** is a 21st Century, State of the Art interface unit with features **ONLY** offered by comPPage. The New SQL Server Database Enterprise system. Adds more functions while maintaining the same easy to use programming interface. Along with the upgrades is a new **Reporting** system for generating custom reports both local and from a remote site. Improvements include additional **Escalations, Building numbers, Wing identification and Floor. Wireless for an easy install.***

System:	Controller, SQL Database, Radmin®, PIA Alerts Display Software, Router.
Connectivity:	Connect PIA units to Controller via RS485, Ethernet, or Wireless.
Accepted Inputs:	Voltage Free, Voltage, and Pulsing or Mixed
Paging Protocol:	Scope, TAP1.8, Comp2 and WaveWare.
System Capability:	AiPhone, Cornell, Engineered Electronics, Heritage MedCall, Tektone Sound & Signal, and Dukane
PIA Approvals:	FCC Title 47 Part 15, ANSI C63.4-2003, UL 1069 Ed: 7, CSA C22.2#205, RoHS Compliant,

Major Features

- ▶ **Ability to Disable Reset Message for Notification but not for logging**
- ▶ **Multiple Paging Outputs – RS232, IP Address, Email and SMS**
- ▶ **Eight Levels of Escalation – set individually for each Zone**
- ▶ **Unlimited number of Paging Interface Adapters attached to Controller**
- ▶ **Wireless Paging Interface Adapters and Router**
- ▶ **Custom Reports – Date/Time, Building, Floor, Wing, Room, Bed Average and Elapse Time Report**



Report Print Out

Alert Notification Report								
Report Date: 01/01/2011			Date: From: 06/29/2010 To: 06/29/2010					
Location: Happy Valley Retirement								
Bldg	Floor	Wing	Room	Bed	Alarm Message	Activated	Reset	Response Time
06/29/2010								
1	3	East	300	A	Bed Call Room 300	03:09:33 PM	03:12:24 PM	00:02:61
1	3	East	300	A	Bath Call Room 300	03:25:49 PM	03:27:43 PM	00:01:54
1	3	East	300	A	Bed Call Room 300	09:31:12 PM	09:31:12 PM	00:08:25
Average								00:04:23
3	1		100		Bath Call Room 100	10:20:44 AM	10:24:07 AM	00:03:23
3	2		200		Bed Call Room 200	02:40:06 PM	02:41:09 PM	00:01:03
3	1		102		Bed Call Room 102	07:57:44 PM	08:28:29 PM	00:30:45
3	1		105		Bed Call Room 105	11:53:06 PM	11:56:59 PM	00:03:53
Average								00:09:46

Service - RS-232 Transmitter on Server - 1

IP Configuration IP Addr: 192.168.1.184 IP Port: 0		Transmitter Settings Port Name: COM1 Port Settings: 9600-N-8-1 Protocol: Scope Data Type: Alphanumeric Baud Rate: 512	
Notifications/Logging Send Notifications: <input checked="" type="checkbox"/> Log Notification Data: <input checked="" type="checkbox"/>			

Controller Maintenance

Connection Edit Add Delete

Unit - 00:90:C2:D3:F4:12

- Zone - 1
- Zone - 2
- Zone - 3
- Zone - 4
- Zone - 5
- Zone - 6
- Zone - 7
- Zone - 8
- Zone - 9
- Zone - 10
- Zone - 11
- Zone - 12
- Zone - 13
- Zone - 14
- Zone - 15
- Zone - 16
- Zone - 17
- Zone - 18
- Zone - 19
- Zone - 20
- Zone - 21
- Zone - 22
- Zone - 23
- Zone - 24
- Zone - 25
- Zone - 26
- Zone - 27
- Zone - 28
- Zone - 29
- Zone - 30

Zone - 1 of Unit - 00:90:C2:D3:F4:12 on Server - 1

Location Wing Number: West Room Number: 100 Bed Number:	Control Info Connection Type: Voltage Alert-High Repeat Interval: 2 Minutes Repeat Count: 3	Overrides Suppress Reset Message: <input type="checkbox"/>																				
Messages Active: Bed Call Room 100 Reset: Room 100 Reset Slow Pulse: <input checked="" type="checkbox"/> Bath Call Room 100 Fast Pulse: <input checked="" type="checkbox"/> Emergency Call Room 100																						
Escalation Settings <table border="1"> <thead> <tr> <th></th> <th>Interval</th> <th></th> <th>Interval</th> </tr> </thead> <tbody> <tr> <td>First: <input checked="" type="checkbox"/></td> <td>10 Minutes</td> <td>Fifth: <input type="checkbox"/></td> <td>2 Minutes</td> </tr> <tr> <td>Second: <input type="checkbox"/></td> <td>15 Minutes</td> <td>Sixth: <input type="checkbox"/></td> <td>2 Minutes</td> </tr> <tr> <td>Third: <input type="checkbox"/></td> <td>2 Minutes</td> <td>Seventh: <input type="checkbox"/></td> <td>2 Minutes</td> </tr> <tr> <td>Fourth: <input type="checkbox"/></td> <td>2 Minutes</td> <td>Eighth: <input type="checkbox"/></td> <td>2 Minutes</td> </tr> </tbody> </table>			Interval		Interval	First: <input checked="" type="checkbox"/>	10 Minutes	Fifth: <input type="checkbox"/>	2 Minutes	Second: <input type="checkbox"/>	15 Minutes	Sixth: <input type="checkbox"/>	2 Minutes	Third: <input type="checkbox"/>	2 Minutes	Seventh: <input type="checkbox"/>	2 Minutes	Fourth: <input type="checkbox"/>	2 Minutes	Eighth: <input type="checkbox"/>	2 Minutes	Notification Settings Pager E-Mail SMS Text
	Interval		Interval																			
First: <input checked="" type="checkbox"/>	10 Minutes	Fifth: <input type="checkbox"/>	2 Minutes																			
Second: <input type="checkbox"/>	15 Minutes	Sixth: <input type="checkbox"/>	2 Minutes																			
Third: <input type="checkbox"/>	2 Minutes	Seventh: <input type="checkbox"/>	2 Minutes																			
Fourth: <input type="checkbox"/>	2 Minutes	Eighth: <input type="checkbox"/>	2 Minutes																			

Notification Settings

Pager

E-Mail

SMS Text

Pager IDs of Zone - 1 of Unit - 00:90:C2:D3:F4:12 on Server - 1

Shifts		Escalations	
Pager IDs:	Beep Types:	Pager IDs:	Beep Types:
Shift1 ID 1: 0300800	Type A	First: 0304000	Type D
Shift1 ID 2: 0301600	Type B	Second: 0000000	(None)
Shift1 ID 3: 0302400	Type C	Third: 0000000	(None)
Shift2 ID 1: 0000000	(None)	Fourth: 0000000	(None)
Shift2 ID 2: 0000000	(None)	Fifth: 0000000	(None)
Shift2 ID 3: 0000000	(None)	Sixth: 0000000	(None)
Shift3 ID 1: 0000000	(None)	Seventh: 0000000	(None)
Shift3 ID 2: 0000000	(None)	Eighth: 0000000	(None)
Shift3 ID 3: 0000000	(None)		