


PROGRAMMING

The Scope® LTX transmitter installed in the comPPage® Staff-Links Bed/Bath Call station has been programmed with a default file. You should **Read** the device and modify it as necessary. Be sure to read the documentation that comes with the unit before Programming and installing the unit.

1. Frequency – 446.525, Baud – 512, Pager Cap Code - 0300800
2. Power
 - a. The power selection is set for the battery shipped with the unit. This should not be Changed. (3Volts DC)
3. Insert the battery observing polarity.
4. If programming a bed station be sure to insert the supplied plastic dummy plug.
5. Install the programming utility software on a computer.
6. Insert the 3V battery.
7. Remove the LED connector from the programming header.
8. Connect the programming cable to a USB Port on the computer and the other end to the unit To be programmed. The end that is attached to the unit has a guide on the side of the header Pin. It should be facing away from the power terminal.
9. After connecting the programming cable bring up the program utility and select [READ]. Once you read the device you can modify the file and [WRITE] it back to the unit. It is Recommended to save your file so you can reload if necessary.
10. Re-install the LED connector to the program header.



Programming Header

- Contact #1 – Bed Call Cord or Bath Pull cord
- Contact #2 – Help Button
- Contact #3 – Staff Assist (under logo)
- Contact #4 – Cancel Button

INSTALLATION

1. Mount the rear plate to the wall using the supplied screws.
2. Carefully slide the unit over the rear plate, taking care not to disrupt the ribbon cable.
3. Bed Station: Remove the plastic dummy plug and insert the pear-push call cord and Push the cancel button. This will cancel the call.
4. Insert the bottom screw that holds the Front part of the unit to the rear plate.



Rear Plate

OPERATION

1. When the **Help button** is pushed the unit will transmit a [Help needed in RM 22].
2. When a **Staff Member** is in the room and needs assistance, there is a hidden button Behind the [c] in comPPage. It is programmed to transmit [Staff Assist in RM22].
3. The Pull String is programmed to transmit [Bath Call RM 22].
4. If a Call Cord is installed it will work the same as the Pull String.
5. All calls are set to repeat till canceled.

