

comPPage<sup>®</sup>, Inc.

PIA Enterprise System

# comPPage<sup>®</sup>, Inc.

Controller w/SQL Database  
PIA and RS232 Input



**Made in The USA**

# SQL Server Software

- \* 64 Contact Inputs per PIA unit
- \* Fused Non Voltage Inputs
- \* Voltage or Non-Voltage and Slow or Fast Pulse Host
- \* AC or DC Input
- \* 5-50 Volts Input
- \* Supports Positive or Negative Common
- \* Eight Levels of Escalation
- \* Supports Shifts (3)
- \* Unlimited Interface Units on one Controller
- \* Supports Different Host systems on the Same Controller
- \* Supports RS232 inbound from Panels and NC systems


# SQL Server Software

- \* Email, SMS Texting
- \* IP & RS232 Pocket Paging
- \* Direct Connect to EnGenius Cordless Telephones
- \* Wireless to Smart Device(Android) on WiFi Local Area Network
- \* Create and Print Custom Reports
- \* Display Alerts on All or Selected monitors on the Local Area Network
  - \* Bed Calls – displayed Black on White
  - \* Bath Calls - displayed Black on **Blue**
  - \* Emergency Calls – displayed Black on **Red**
  - \* Call Escalations – displayed Black on **Green**
- \* Display Alerts on HD TV (Option)

# PIA3264R2

- \* UL 1069 Standard for Safety Hospital Signaling and Nurse Call Equipment Issue:2007/10/12Ed:7
- \* CSA C22.2#205 Signaling Equipment General Instruction No 1. Issue: 1983/01/01 (R2004)
- \* FCC – CFR 47, Part 15
- \* RoHS Compliant

# PIA3264R2



- iPhone
- Cornell
- Dukane
- Engineered Electronics
- Heritage MedCall
- Jeron
- Tektone Sound + Signal

# Enterprise System

**Wireless Operation is an option**

- iPhone
- Cornell
- Heritage Medcall

- Engineered Electronics
- Tektone Sound + Signal



- 457 Mhz
- Btn Alert
  - Cord Alert
  - Fall Alert
  - Test Call
  - Low battery



Relay Control



Wireless Router


- Internet Connection
- SMS Texting
- Emailing



IP Paging



EnGenius  
Telephones  
TAP 1.8 Input



Controller

Remote Monitoring of Alerts on any hardwired or wireless Network Station



### Alert Notification Report


Custom Reports

Report Date: 3/28/2011  
 Location: Happy Valley Retirement

Bldg	Wing	Floor	Room	Bed	Alarm Message	Activated	Reset	Response Time
3/25/2011								
SOUTH	EAST	6	500	A	Bed Call Room 500A	5:53:51 AM	05:54:07 AM	00:00:16
SOUTH	EAST	6	500	A	Bed Call Room 500A	7:17:02 AM	07:17:15 AM	00:00:13
SOUTH	EAST	6	500	A	Bath Call Room 500	7:17:24 AM	07:17:33 AM	00:00:09
SOUTH	EAST	6	500	A	Bed Call Room 500A	9:03:23 AM	09:03:29 AM	00:00:06
Total Calls:						4	Average: 00:00:11	

# Enterprise System

## SERIAL NETWORK ALERT PROCESSOR





# Digital Alert Viewer Paging

## HDTV



- ▶ **Display Alerts on HD TV**
- ▶ **Option of 9 different Colors**
- ▶ 7 Cap Codes Available
- ▶ Use on any POCSAG Transmitter
- ▶ Frequencies – 456-459, 464-467
- ▶ Connect via HDMI cable

